

PURIKER RK

Experimental rocker with durable and smooth speed change




STANDARD SPECIFICATION

MODEL	PURIKER RK
Speed range	4 ~ 80 rpm
Tilt Angle	±9°
Motor	BLDC Motor
Speed Control	Analog Jog Dial
Platform size (WxDxH)mm	300 x 300
Size (WxDxH)mm	300 x 378 x 142
Power Consumption	20W
Power	AC100 ~ 240V, 50/60Hz
Weight(Kg)	5.4
Cat. No.	CS-PMX-Set

MAIN FEATURE

- Seesaw motion [2D motion]
- Equipped with a long-life BLDC motor without noise, spark, and electrical noise due to friction
- Silicone pad as standard
2nd level platform(option)



ORDER INFORMATION

DIVISION	ITEM	PURIKER RK
		Cat. No.
PRODUCT	Main Set	CS-PKR-RK-Set